

Recombinant Human Kallikrein 4, KLK4 (C-6His)

Catalog No: #AP75083



Package Size: #AP75083-1 10ug #AP75083-2 50ug #AP75083-3 500ug #AP75083-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Kallikrein 4, KLK4 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9Y5K2
Uniprot	Q9Y5K2
GeneID	9622;
Calculated MW	25.44kDa
Target Sequence	SCSQIINGEDCSPHSQPWQAALVMENELFCSGVLVHPQWVLSAAHCFQNSYTIGLGLHSLEADQEPGSQMVE ASLSVRHPEYNRPLLANDLMLIKLDESSESMTIRSISIASQCPTAGNSCLVSGWGLLANGRMPTVLQCVNVS VSEEVCSKLYDPLYHPSMFCAGGGQDQKDCSNGDGGPLICNGYLQGLVSGKAPCGQVGVPGVYTNLCKF TEWIEKTVQASVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only