

Recombinant Human Kallikrein 8,Neuropsin (C-6His)

Catalog No: #AP75087



Package Size: #AP75087-1 10ug #AP75087-2 50ug #AP75087-3 500ug #AP75087-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Kallikrein 8,Neuropsin (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O60259
Uniprot	O60259
GeneID	11202;
Calculated MW	26.04kDa
Target Sequence	QEDKVLGGHECQPHSQPWQAALFQGQQLLCGGVLVGGNWWLTAAHCKKPKYTVRLGDHSLQNKDGPQEI PVVQSIPHPCYNSSDVEDHNHDLMLLQLRDQASLGSVKPISLADHCTQPGQKCTVSGWGTVTSPRENF PDT LNCAEVKIFPQKKCEDAYPGQITDVMVCAGSSKGADTCQGDSSGGLVCDGALQGITSWGS DPCGRSDKPGV YTNICRYLDWIKKIIGSKGVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only