

Recombinant Human Legumain,Asparaginyl Endopeptidase (C-6His)



Catalog No: #AP75091

Orders: order@signalwayantibody.com

Package Size: #AP75091-1 10ug #AP75091-2 50ug #AP75091-3 500ug #AP75091-4 1mg

Support: tech@signalwayantibody.com

Description

| | |
|-----------------|---|
| Product Name | Recombinant Human Legumain,Asparaginyl Endopeptidase (C-6His) |
| Host Species | Human Cells |
| Purification | Greater than 95% as determined by reducing SDS-PAGE. |
| Specificity | Human |
| Accession No. | Q99538 |
| Uniprot | Q99538 |
| GeneID | 5641; |
| Calculated MW | 48.72kDa |
| Target Sequence | IPIDDPEDGGKHWVIVAGSNGWYNRYHQADACHAYQIIHRNGIPDEQIVMMYDDIAYSEDNPTPGIVINRPN GTDVYQGVPKDYTGEDVTPQNFLAVLRGDAEAVKIGISGKVLKSGPQDHFVFIYFTDHGSGTILVFPNEDLHVK DLNETIHYMYKHKMYRKMVFYIEACESGSMNHLDPNINVYATTAANPRESSYACYDEKRSTYLGDWYSVN WMEDSDVEDLTKETLHKQYHLVKSHTNTSHVMQYGNKTISTMKVMQFQGMKRKASSPVPLPPVTHLDLTPS PDVPLTIMKRKLMNTNDLEESRQLTEEIQRHLDARHLIEKSVRKIVSLLAASEAEVEQLLSERAPLTGHSCYPEA LLHFRTHCFNWHSPTYEYALRHLYVLVNLCEKPYPLHRIKLSMDHVCLGHYVDHHHHHH |
| Formulation | Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 10% Glycerol, pH 8.0. |
| Storage | Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles. |

Note: This product is for in vitro research use only