

Recombinant Human Matrix Metalloproteinase-1, MMP-1 (C-6His)

Catalog No: #AP75093

Orders: order@signalwayantibody.com

Package Size: #AP75093-1 10ug #AP75093-2 50ug #AP75093-3 500ug #AP75093-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Matrix Metalloproteinase-1, MMP-1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P03956
Uniprot	P03956
GeneID	4312;
Calculated MW	52.88kDa
Target Sequence	FPATLETQEVDVLDVQKYLEKYNNLKNDRQVEKRRNSGPVVEKLGKQMQEFFGLKVTGKPD AETLKV MKQP RCGVDPVAQFVLTEGNPRWEQTHLTYRIENYTPDLPRADV DHAIEKAFQLWSNV TPLFTTKVSEGGADIMISF VRGDHRDNSPFDGPGGNLAHAFQPGPGIGGDAHFDEDERWTN NFREYNLHRVAAHELGHSLGLSHSTDIGA LMYPSYTFSGDVQLAQDDIDGIAIYGRSQNPVQPIGPQTPKACDSKLTFDAIT TIRGEVMFFKDRFYMR TNPF YPEVELNFISVFWPQLPNGLEAAYEFADRDEV RFFKGNKYWAVQGGQNVLHGYPKDIYSSFGFPRTVKHIDAAL SEENTGKTYFFVANKYWR YDEYKRSM DPGYPK MIAHDFPGIGHKVD AVFMKDGFFYFFH GTRQYKFDPKTKR ILTLQKANSWFNCRKNVDH HHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM MES, 150mM NaCl, 0.05% Brij35, pH 5.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only