

# Recombinant Human Matrix Metalloproteinase-3, MMP-3 (C-6His)



Catalog No: #AP75095

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75095-1 10ug #AP75095-2 50ug #AP75095-3 500ug #AP75095-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Matrix Metalloproteinase-3, MMP-3 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P08254
Uniprot	P08254
GeneID	4314;
Calculated MW	53.26kDa
Target Sequence	YPLDGAARGEDTSMNLVQKYLENYDLEKDVKQFVRRKDSGPVVKKIREMQKFLGLEVTGKLDSDTLEVMRK PRCGVPDVGHFRTFPGIPKWRKTHLTYRIVNYTPDLPKDAVDSAVEKALKVWEEVTPLTFSRLYEGEADIMISF AVREHGDFYPFDGPGNVLAHAYAPGPGINGDAHFDDEQWTKDTTGTNLFLVAAHEIGHSLGLFHSANTEAL MYPLYHSLDLTRFRLSQDDINGIQSLYGPPDSETPPLVPTPEPVPPEPGTPANCDPALSFDVAVSTLRGEILIFK DRHFWRKSLRKLEPELHLISSFWPSLPSGVDAAYEVTSKDLVFIFKGNQFWAIRGNEVRAGYPRGIHTLGFPF TVRKIDAAISDKEKNKTYFFVEDKYWRFDEKRNSMEPGFPKQIAEDFPGIDSKIDAVFEEFGFFYFFTGSSQLE FDPNAKKVTHTLKSNSWLNCDVHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 0.05% Brij35, 10% Glycerol, pH 7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only