

Recombinant Human Polymeric Immunoglobulin Receptor, PIgR (C-6His)



Catalog No: #AP75098

Orders: order@signalwayantibody.com

Package Size: #AP75098-1 10ug #AP75098-2 50ug #AP75098-3 500ug #AP75098-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Polymeric Immunoglobulin Receptor, PIgR (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P01833
Uniprot	P01833
GeneID	5284;
Calculated MW	68.88kDa
Target Sequence	KSPIFGPEEVNSVEGNSVSITCYPPTSVMNRHTRKYWCRQGARGGCITLISSEGYVSSKYAGRANLTNFPENG TFVVNIAQLSQDDSGRYKCGLGINSRGLSFDVSLEVSQGPGLLNDTKVYTVDLGRTVTINCPFKTENAKRKS LYKQIGLYPVLVIDSSGYVNPNYTGRIRLDIQGTGQLLFSVVINQLRLSDAGQYLCQAGDSDNSNKKNADLQVL KPEPELVYEDLRGSVTFHCALGPEVANVAKFLCRQSSGENCDVVVNTL GKRAPAFEGRILLNPQDKDGSFSV VITGLRKEDAGRYLCGAHS DGQLQEGSPIQAWQLFVNEESTIPRSPTVVKGVAGGSVAVLCPYNRKESKSIKY WCLWEGAQNGRCPLLVDSEGWVKAQYEGRLSLEEPGNGTFTVILNQLTSRDAGFYWCLTNGDTLWRTTVE IKIIEGEPNLKVPGNVTAVLGETLKVPCHPCKFSSYEKYWCKWNNTGCQALPSQDEGPSKAFVNCDENSRLV SLTLNLVTRADEGWYWCGVKQGHFYGETAAVYVAVEERKAAGSRDVS LAKADAAPDEKVLDSGFREIENKAI QDPRLFAEEKAVADTRDQADGSRASVDSGSSEEQGGSSRVDH HHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only