

Recombinant Human SLAM Family Member 6,SLAMF6,CD352,NTB-A (C-6His)

Catalog No: #AP75102

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75102-1 10ug #AP75102-2 50ug #AP75102-3 500ug #AP75102-4 1mg

Description

Product Name	Recombinant Human SLAM Family Member 6,SLAMF6,CD352,NTB-A (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q96DU3
Uniprot	Q96DU3
GeneID	114836;
Calculated MW	23.98kDa
Target Sequence	LMVNGILGESVTLPLEFPAGEKVNFIWLFNETSLAFIVPHETKSPEIHVTNPKQGKRLNFTQSYSLQLSNLKME DTGSYRAQISTKTSAKLSSYTLRILRQLRNIQVTNHSQLFQNMTCELHLTCSVEDADDNVSFRWEALGNTLSS QPNLTVSWDPRISSEQDYTCIAENAVSNLSFSVSAQKLCEDVKIQYTDTKVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only