

Recombinant Human Cathepsin B,CTSB (C-6His)

Catalog No: #AP75113



Package Size: #AP75113-1 10ug #AP75113-2 50ug #AP75113-3 500ug #AP75113-4 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Cathepsin B,CTSB (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P07858
Uniprot	P07858
GeneID	1508;
Calculated MW	36.9kDa
Target Sequence	RSRPSFHPVSDELVNYVNKRNTTWQAGHNFYNVDMGYLKRLCGTFLGGPKPPQRMFTEDLKLPASFDARE QWPQCPTIKEIRDQGSCGSCWAFGAVEAISDRICIHTNAHVSVEVSAEDLLTCCGSMCGDGCNGGYPAEAWN FWTRKGLVSGGLYESHVGCPRYSIPPCEHHVNGSRPPCTGEGDTPKCSKICEPGYSPTYKQDKHYGYNYS VSNSEKDIMAEIYKNGPVEGAFSVYSDFLLYKSGVYQHVTGEMMGGHAIRILGWGVENGTPYWLVANSWNTD WGDNGFFKILRGQDHCGIESEVVAGIPRTDQYWEKIVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only