

Recombinant Human Cathepsin E,CTSE (C-6His)

Catalog No: #AP75114



Package Size: #AP75114-1 10ug #AP75114-2 50ug #AP75114-3 500ug #AP75114-4 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Cathepsin E,CTSE (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P14091
Uniprot	P14091
GeneID	1510;
Calculated MW	41.78kDa
Target Sequence	SLHRVPLRRHPSLKKKLRARSQLSEFWKSHNLDMIQFTESCSDQSAKEPLINYLDMEYFGTISIGSPPQNFTV IFDTGSSNLWVPSVYCTSPACKTHSRFQPSQSSTYSQPGQSFSIQYGTGSLSGIIGADQVSVEGLTVVGQQFG ESVTEPGQTFVDAEFDGILGLGYPSLAVGGVTPVFDNMMAQNLVDLPMFSVYMSSNPEGGAGSELIFGGYDH SHFSGSLNWVPVTKQAYWQIALDNIQVGGTVMFCSEGCQAIVDTGTSLITGPSDKIKQLQNAIGAAPVDGEYA VECANLNVMPDVTFTINGVPYTLSPYATLLDFVDGMQFCSSGFQGLDIHPPAGPLWILGDVFIRQFYSVFD RG NNRVGLAPAVPVDHHHHHH
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM MES, 150mM NaCl, pH 5.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only