

Recombinant Human Glypican 3, GPC3, OCI5 (C-6His)

Catalog No: #AP75120



Package Size: #AP75120-1 10ug #AP75120-2 50ug #AP75120-3 500ug #AP75120-4 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Glypican 3, GPC3, OCI5 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P51654
Uniprot	P51654
GeneID	2719;
Calculated MW	61.6kDa
Target Sequence	QPPPPPPDATCHQVRSFFQRLQPGLKWVPETPVPGSDLQVCLPKGPTCCSRKMEEKYQLTARLNMEQLLQS ASMELKFLIIQNAAVFQEAFEIVVRHAKNYTNAMFKNNYPSLTPQAFEFVGEFFTDVSLYILGSDINVDDMVNEL FDSLFPVIYTQLMNPGLPDSALDINECLRGARRDLKVFGNFPKLIMTQVSKSLQVTRIFLQALNLGIEVINTDHL KFSKDCGRMLTRMWYCSYCQGLMMVKPCGGYCNVVMQGCMAGVVEIDKYWREYILSLEELVNGMYRIYDM ENVLLGLFSTIHDSIQYVQKNAGKLTTTIGKLCAHSQQRQYRSAYYPEDLFIDKKVLKVAHVEHEETLSSRRREL IQKLKSFISFYALPGYICSHSPAENDTLCWNGQELVERYSQKAARNGMKNQFNLHELKMGPEPVVSQIIDK LKHINQLLRTMSMPKGRVLDKNLDEEGFESGDCGDEDECIGSGDGMIVKNQLRFLAELAYDLVDVDDAPG NSQQATPKDNEISTFHNLGNVHHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only