

# Recombinant Human Prostatic Acid Phosphatase,ACPP (C-6His)



Catalog No: #AP75125

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75125-1 10ug #AP75125-2 50ug #AP75125-3 500ug #AP75125-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Prostatic Acid Phosphatase,ACPP (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P15309
Uniprot	P15309
GeneID	55;
Calculated MW	42.01kDa
Target Sequence	KELKFVTLVFRHGDRSPIDTFPTDPIKESSWPQGFQQLTLGMEQHYELGEYIRKRYRKFLNESYKHEQVYIR STDVDRTLMSAMTNLAALFPPEGVSIWNPILLWQPIPVHTVPLSEDQLLYLPFRNCPRFQELESETLKSEEFQK RLHPYKDFIATLGKLSGLHGQDLFGIWSKVYDPLYCESVHNFTLPSWATEDMTKLRELSLSSLYGIHKQK EKSRLQGGVLVNEILNHMKRATQIPSYKKLIMYSAHDTTVSGLQMALDVYNGLLPPYASCHLTELYFEKGEYFV EMHYRNETQHEPYPLMLPGCSPSCPLERFAELVGPVIPQDWSTECMTTNSHQGTEDSTDVDHSHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only