

Recombinant Human Acrosomal Protein
SP-10,ACRV1 (C-6His)

Catalog No: #AP75126



Package Size: #AP75126-1 10ug #AP75126-2 50ug #AP75126-3 500ug #AP75126-4 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Acrosomal Protein SP-10,ACRV1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P26436
Uniprot	P26436
GeneID	56;
Calculated MW	26.9kDa
Target Sequence	QPNELSGSIDHQTSVQQLPGEFFSLENPSDAEALYETSSGLNTLSEHGSSEHGSSKHTVAEHTSGEHAESEH ASGEPAATEHAEGEHTVGEQPSGEHLSGEQPLSELESGEQPSDEQPSGEHGSGEQPSGEQASGEQ PSGEHASGEQASGAPISSTSTGTILNCYTCAYMNDQGKCLRGEGETCITQNSQQCMLKKIFEKGKLQFMVQGC ENMCPSMNLFSHGTRMQIICCRNQSFCKIVDHHHHHH
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only