

Recombinant Human Anterior Gradient Protein 2
Homolog,AG-2,HPC8,AGR2 (C-6His)



Catalog No: #AP75128

Package Size: #AP75128-1 10ug #AP75128-2 50ug #AP75128-3 500ug #AP75128-4 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Anterior Gradient Protein 2 Homolog,AG-2,HPC8,AGR2 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O95994
Uniprot	O95994
GeneID	10551;
Calculated MW	18.85kDa
Target Sequence	RDTTVKPGAKKDTKDSRPKLPQTLSRGWGDQLIWTQTYEEALYKSKTSNKPLMIIHHLDECPHSQALKKVFAE NKEIQKLAEQFVLLNLVYETTDKHLSPDGQYVPRIMFVDPSLTVRADITGRYSNRLYAYEPADTALLLDNMKKA LKLLKTELVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris,200mM NaCl,10GlycerolpH8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only