

Recombinant Human
Fetuin-A,AHSG, α -2-HS-Glycoprotein,AHSG (C-6His)

Catalog No: #AP75129

Package Size: #AP75129-1 10ug #AP75129-2 50ug #AP75129-3 500ug #AP75129-4 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Fetuin-A,AHSG, α -2-HS-Glycoprotein,AHSG (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P02765
Uniprot	P02765
GeneID	197;
Calculated MW	38.36kDa
Target Sequence	APHGPGLIYRQPNCCDPETEEAALVAIDYINQNLPGYKHTLNQIDEVKVWPQQPSGELFEIEIDTLETTCHVL DPTPVARCSVRQLKEHAVEGDCDFQLLKLDGKFSVVYAKCDSSPDSAEDVRKVCQDCPLLAPLNDTRVVHAA KAALAAFNAQNNGSNFQLEEISRAQLVPLPPSTYVEFTVSGTDCVAKEATEAAKCNLLAEKQYGFCKATLSEK LGGAEVAVTCTVFQTQPVTSQPQPEGANEAVTPVVDPDAPSPPLGAPGLPPAGSPPDShVLLAAPPQHQL HRAHYDLRHTFMGVVSLGSPSGEVSHPRKTRTVVQPSVGAAAGPVVPPCPGRIRHFKVVDHHHHHH
Formulation	Lyophilized from a 0.2 μ m filtered solution of 20mM TrisHCl, 150mM NaCl, pH 7.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only