



Recombinant Human Asialoglycoprotein Receptor 1,ASGPR1 (C-6His)

Catalog No: #AP75131

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #AP75131-1 10ug #AP75131-2 50ug #AP75131-3 500ug #AP75131-4 1mg

Description

Product Name	Recombinant Human Asialoglycoprotein Receptor 1,ASGPR1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P07306
Uniprot	P07306
GeneID	432;
Calculated MW	27.43kDa
Target Sequence	QNSQLQEELRGLRETFSNFTASTEAQVKGLSTQGGNVGRKMKSLESQLEKQQKDLSEDHSSLLLVHKQFVSD LRSLSCQMAALQGNGSERTCCPVNWVEHERSCYWFSRSGKAWADADNYCRLEDAHLVVVTSWEEQKFVQ HHIGPVNTWMGLHDQNGPWKWVDGTDYETGFKNWRPEQPDDWYGHGLGGGEDCAHFTDDGRWNDDVC QRPYRWVCETELDKASQEPPLLVDHHHHHH
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only