

Recombinant Human Brain-Specific Angiogenesis Inhibitor 3,BAI3 (C-6His)



Catalog No: #AP75134

Package Size: #AP75134-1 10ug #AP75134-2 50ug #AP75134-3 500ug #AP75134-4 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Brain-Specific Angiogenesis Inhibitor 3,BAI3 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O60242
Uniprot	O60242
GeneID	577;
Calculated MW	97.51kDa
Target Sequence	AQDFWCSTLVKGVYGSYSVSEMFPKNFTNCTWTLENPDPTKYSIYLKFSKKDLSCSNFSLLAYQFDHFSHEKI KDLLRKNHSIMQLCNSKNAFVFLQYDKNFIQIRRVFPTNFPGLQKKGEEDQKSFFEFLVLNKVSPSQFGCHVLC TWLESLCKSENGRTEGIMYTKCTCPQHLGEWGIDDQSLILLNNVVLPLNEQTEGCLTQELQTTQVCNLTRE AKRPPKEEFGMMGDHTIKSQRPVHEKRVPEQADAAKFMAQTGESGVEEWSQWSTCSVTGCGGSQVR TRTCVSPYGTHCSGPLRESRVCNNTALCPVHGVWEEWSPWSLCSFTCGRGQRTRTRRSCTPPQYGGRPCEG PETHHKPCNIALCPVDGQWQEWSSWSQCSVTCSNGTQQRSRQCTAAAHGGSECRGPWAESRECYNPECT ANGQWNQWGHWSGCSKSCDGGWERRIRTCQGAVITGQCEGTGEEVRRCSERCPAPYEICPEDYLMMSM VWKRTAGDLAFNQCPLNATGTTSRRCSLSLHGVAFWEQPSFARCISNEYRHLQHSIKEHLAKGQRMLAGDG MSQVTKTLLDLTQRKNFYAGDLLMSVEILRNVTDTFKRASYIPASDGVQNFFQIVSNLLDEENKEKWEDAQQIY PGSIELMQVIEDFIHIVGMGMMDFQNSYLMTGNVVASIQKLPAAVSLTDINFPMKGRKGMVDWARNSEDRVVI PKSIFTVPVSSKELDESSVFLGAVLYKNLDLILPTLRNYTVINSKIIVTIRPEPKTTDSFLEIELAHLANGTLNPYC VLWDDSKTNESLGTWSTQGCKTVLTDASHTKCLCDRLSTFAILAAQPREIIMESSGTPSVTVDHSHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 5% Thehalose, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only