

Recombinant Human EMMPRIN,Basigin,CD147
(C-6His)

Catalog No: #AP75136

Package Size: #AP75136-1 10ug #AP75136-2 50ug #AP75136-3 500ug #AP75136-4 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human EMMPRIN,Basigin,CD147 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P35613-2
Uniprot	P35613
GeneID	682;
Calculated MW	21.16kDa
Target Sequence	AAGTVFTTVEDLGSKILLTCSLNSATEVTGHRWLKGGVVLKEDALPGQKTEFKVDSDDQWGEYSCVFLPEP MGTANIQLHGPPRVKAVKSSEHINEGETAMLVCKSESVPVTDWAWYKITDSEDKALMNGSESRRFFVSSSQG RSELHIENLNMEADPGQYRCNGTSSKGSDQAIITLRVRSHVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only