

Recombinant Human Carnosine Dipeptidase 1,CNDP1 (C-6His)

Catalog No: #AP75139

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75139-1 10ug #AP75139-2 50ug #AP75139-3 500ug #AP75139-4 1mg

Description

Product Name	Recombinant Human Carnosine Dipeptidase 1,CNDP1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q96KN2
Uniprot	Q96KN2
GeneID	84735;
Calculated MW	54.9kDa
Target Sequence	<p>PSPPPALLEKVFQYIDLHQDEFVQTLKEWVAIESDSVQPVPRFRQELFRMMAVAADTLQRLGARVASVDMGP</p> <p>QQLPDGQSLPIPPVILAEELGSDPTKGTVCFYGHLDVQPADRGDGLWLTDPYVLTEVDGKLYGRGATDNKGPVL</p> <p>AWINAVSAFRALEQDLPVNIKFIIIEGMEEAGSVALEELVEKEKDRFFSGVDYIVISDNLWISQRKPAITYGTRGNS</p> <p>YFMVEVKCRDQDFHSGTFGGILHEPMADLVALLGSLVDSSGHILVPGIYDEVVPLTEEEINTYKAIHLDLEEYRN</p> <p>SSRVEKFLFDTKEEILMHLWRYPSLSIHGIEGAFDEPGTKVIPGRVIGKFSIRLVPHMNVSAVEKQVTRHLEDV</p> <p>FSKRNSSNMVVSMTLGLHPWIANIDDTQYLAAKRAIRTVFGTEPDMIRDGSTIPIAKMFQEIVHKSVVLIPLGA</p> <p>VDDGEHSQNEKINRWNYIEGTKLFAAFFLEMAQLHVDHHHHHH</p>
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 7.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only