



Recombinant Human Fc γ RIIIA,FCGR3A,CD16a (C-6His)

Catalog No: #AP75143

Package Size: #AP75143-1 10ug #AP75143-2 50ug #AP75143-3 500ug #AP75143-4 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Fc γ RIIIA,FCGR3A,CD16a (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P08637
Uniprot	P08637
GeneID	2214;
Calculated MW	22.61kDa
Target Sequence	GMRTEDLPKAVVFLEPQWYRVLEKDSVTLKCQGAYSPEDNSTQWFHNESLISSQASSYFIDAATVDDSGEYR CQTNLSTLSDPVQLEVHIGWLLLQAPRWVFKEEDPIHLRCHSWKNTALHKVTYLQNGKGRKYFHHNSDFYIPK ATLKDSGSYFCRGLVGSKNVSSETVNITITQGLAVSTISSFFPPGYQH HHHHHH
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only