

Recombinant Human CD99, MIC2 (C-6His)

Catalog No: #AP75149



Package Size: #AP75149-1 10ug #AP75149-2 50ug #AP75149-3 500ug #AP75149-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CD99, MIC2 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P14209
Uniprot	P14209
GeneID	4267;
Calculated MW	11.13kDa
Target Sequence	DGGFDLSDALPDNENKKPTAIPKKPSAGDDFDLGDVVDGENDDP RPPNPPKMPNPNPNHPSSSGSFSDA DLADGVSGGEGKGGSDGGGSHRKEGEEADV DHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only