

Recombinant Human Carboxylesterase 1,CES1 (C-6His)

Catalog No: #AP75152

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75152-1 10ug #AP75152-2 50ug #AP75152-3 500ug #AP75152-4 1mg

Description

Product Name	Recombinant Human Carboxylesterase 1,CES1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P23141-3
Uniprot	P23141
GeneID	1066;
Calculated MW	61.05kDa
Target Sequence	HPSSPPVVDTVHGKVLGKFVSLEGFAQPVAIFLGIPFAKPPLGPLRFTPPQPAEPWSFVKNATSYPPMCTQDP KAGQLLSELFTRKENIPLKLSCLYLNITYTPADLTKKNRLPVMVWIHGGGLMVGAASTYDGLALAAHENVV VTIQYRLGIWGFSTGDEHSRGNWGHLDQVAALRWVQDNIAFSGGNPGSVTIFGESAGGESVSVLVLSPLAK NLFHRAISESGVALTSVLVKKGDVKPLAEQIAITAGCKTTTSAVMVHCLRQKTEEELETTLKMKFLSLDLQGD RESQPLLGTVIDGMLLLKTPEELQAERNFHTVPYMGVINKQEFGLIPMLMSYPLSEGQLDQKTAMSLWWSY PLVCIKELIPEATEKYLGGTDDTVKKKDLFLDLIADVMFGVPSVIVARNHRDAGAPTYMYEFQYRPSFSSDMK PKTVIGDHGDELFSVFGAPFLKEGASEEEIRLSKMVMKFWANFARNGNPNGEGLPHWPEYNQKEGYLQIGAN TQAAQKLKDEVAFWTNLFAKKAVEKPPQTEVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM HAc-NaAc, 150mM NaCl, pH 4.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only