

# Recombinant Human Carboxypeptidase B,CPB1 (C-6His)

Catalog No: #AP75157

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75157-1 10ug #AP75157-2 50ug #AP75157-3 500ug #AP75157-4 1mg

## Description

Product Name	Recombinant Human Carboxypeptidase B,CPB1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P15086
Uniprot	P15086
GeneID	1360;
Calculated MW	46.94kDa
Target Sequence	HHGGEHFEGEKVFRVNVEDENHINIIRELASTTQIDFWKPDSVTQIKPHSTVDFRVKAEDTVTVENVLKQNELQ YKVLISNLRNVVEAQFDSRVRATGHSYEKYNKWETIEAWTQQVATENPALISRSVIGTTFEGRAIYLLKVGKAG QNKPAIFMDCGFHAREWISPAFCQWVREAVRTYGREIQVTELLDKLDFYVLPVLPVLDGYYITWTKSRFWRKT RSTHTGSSCIGTDPNRNFDAGWCEIGASRNPCDETYCGPAAESEKETKALADFIRNKLSSIKAYLTIHSYSQMM IYPYSYAYKLGENNAELNALAKATVKELASLHGTYGPGATTIYPAAGGSDWAYDQGIKIRYSFTFELRDTG RYGFLLPESQIRATCEETFLAIKYVASYVLEHLYVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM MES, 150mM NaCl, pH 5.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only