

Recombinant Human Cathepsin L2,CTSL2 (C-6His)

Catalog No: #AP75162



Package Size: #AP75162-1 10ug #AP75162-2 50ug #AP75162-3 500ug #AP75162-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Cathepsin L2,CTSL2 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O60911
Uniprot	O60911
GeneID	1515;
Calculated MW	36.69kDa
Target Sequence	VPKFDQNLDTKWYQWKATHRRLYGANEEGWRRRAVWEKNMKMIELHNGEYSQGKHGFTMAMNAFGDMTNE EFRQMMGCFRNQKFRKGGKVFREPLFLDLPKSVDRKKGYYVTPVKNQKQCGSCWAFSATGALEGQMFRKTG KLVSLSEQNLVDCSRPQGNQGCNGGFMARAFQYVKENGGLDSEESYPYVAVDEICKYRPENSVANDTGFTV VAPGKEKALMKAVATVGPISVAMDAGHSSFQFYKSGIYFEPDCSSKNDHGVLVVGYGFEGANSNNSKYWLV KNSWGPEWGSNGYVKIAKDKNNHCGIATAASYPNVVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM MES, 150mM NaCl, pH 5.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only