

Recombinant Human β -Galactosidase, GLB1 (C-6His)

Catalog No: #AP75175



Package Size: #AP75175-1 10ug #AP75175-2 50ug #AP75175-3 500ug #AP75175-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human β -Galactosidase, GLB1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P16278
Uniprot	P16278
GeneID	2720;
Calculated MW	74.63kDa
Target Sequence	LRNATQRMFEIDYSRDSFLKDGQPFRIYISGSIHYSRVPRFYWKDRLLKMKMAGLNAIQTYVPWNFHEPWPGQ YQFSEDHDVEYFLRLAHELGLLVILRPGPYICAEWEMGGLPAWLLEKESILLRSSDPDYLAAVDKWLGVLLPKM KPLLYQNGGPVITVQVENEYGSYFACDFDYLRFQKRFRHHLGDDVVLFTTDGAHKTFKCGALQGLYTTVDF GTGSNITDAFLSQRKCEPKGPLINSEFYTGWLDHWGQPHSTIKTEAVASSLYDILARGASVNLVYFVIGGTFNFAY WNGANSPYAAQPTSVDYDAPLSEAGDLTEKYFALRNIIQKFEKVPEGPIPPSTPKFAYGKVTLEKLTVGAALDI LCPSGPIKSLYPLTFIQVKQHYGFVLYRRTLPQDCSNPAPLSSPLNGVHDRAYVAVDGIPQGVLERNNVITLNI GKAGATDLLVENMGRVNYGAYINDFKGLVSNLTLSSNILDWTIFPLDTEDAVRSHLGGWGHDRDSGHHDEA WAHSSNYTLPAFYMGNFISIPSGIPDLPQDTFIQFPGWTKGQVWINGFNLGRYWPARGPQLTLFVPQHILMST APNTITVLELEWAPCSSDDPELCAVTFVDRPVISSVTYDHPSPKPVKRLMPPPPQKNKDSWLDHVVDHHHH HH
Formulation	Supplied as a 0.2 μ m filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only