

# Recombinant Human Glycoprotein A33,GPA33,GPA33 (C-6His)



Catalog No: #AP75176

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75176-1 10ug #AP75176-2 50ug #AP75176-3 500ug #AP75176-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Glycoprotein A33,GPA33,GPA33 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q99795
Uniprot	Q99795
GeneID	10223;
Calculated MW	24.66kDa
Target Sequence	ISVETPQDVLRLASQGKSVTLPLCTYHTSTSSREGLIQWDKLLLTHTERVVIWPFNSKNYIHGELYKNRVSISNNA EQSDASITIDQLTMADNGTYECSVSLMSDLEGNTKSRVLLLVPPSKPECGIEGETIIGNNIQLTCQSKEGSPT PQYSWKRYNILNQEQLAQPASGQPVSLKNISTDTSGYYICTSSNEEGTQFCNITVAVRSPSMNVVDHHHHH H
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only