

## Recombinant Human B7-H2,ICOSLG,CD275 (C-6His)

Catalog No: #AP75177



Package Size: #AP75177-1 10ug #AP75177-2 50ug #AP75177-3 500ug #AP75177-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human B7-H2,ICOSLG,CD275 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O75144
Uniprot	O75144
GeneID	23308;
Calculated MW	27.72kDa
Target Sequence	DTQEKEVRAMVGSDELSCACPEGSRFDLNDVYVYWQTSESKTVVYHIPQNSSLENVDSRYRNRALMSPA GMLRGDFSLRLFNVTPQDEQKFHCLVLSQSLGFQEVLSVEVTLHVAANFVSVVVSAPHSPSQDELFTCTCTIN GYPRPNVYWINKTDNSLLDQALQNDTVFLNMRGLYDVVSVLRIARTPSVNI GCCIENVLLQQNLTVGSQTGNDI GERDKITENPVSTGEKNAATWSVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only