

# Recombinant Human Interferon $\alpha$ , $\beta$ Receptor 2,IFNAR2 (C-6His)

Catalog No: #AP75178

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75178-1 10ug #AP75178-2 50ug #AP75178-3 500ug #AP75178-4 1mg

## Description

Product Name	Recombinant Human Interferon $\alpha$ , $\beta$ Receptor 2,IFNAR2 (C-6His)
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P48551
Uniprot	P48551
GeneID	3455;
Calculated MW	25.79kDa
Target Sequence	ISYDSPDYTDESCCTFKISLRNFRSILSWELKNHSIVPTHYLLYTIMSKPEDLKVVKNCANTTRSFCDLTDEWRS THEAYVTVLEGFSGNTTLFSCSHNFWLAIDMSFEPPEFEIVGFTNHINVMVKFPSIVEEELQFDLSLVIEEQSEGI VKKHKPEIKGNMSGNFTYIIDKLIPNTNYCVSVYLEHSDEQAVIKSPLKCTLLPPGQESSESAAESAKVDHHHHHH
Formulation	Lyophilized from a 0.2 $\mu$ M filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only