

Recombinant Human Kallikrein 1, KLK1 (C-6His)

Catalog No: #AP75183



Package Size: #AP75183-1 10ug #AP75183-2 50ug #AP75183-3 500ug #AP75183-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Kallikrein 1, KLK1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P06870
Uniprot	P06870
GeneID	3816;
Calculated MW	28.15kDa
Target Sequence	PPIQSRIVGGWECEQHSQPWQAALYHFSTFQCGGILVHRQWVLTAAHCISDNYQLWLGRHNLFDDENTAQFV HVSESFPHPGFNMSLLENHTRQADEDYSHDLMLLRLTEPADTITDAVKVVELPTQEPEVGSTCLASGWGSI ENFSFPDDLQCVDLKILPNDECKKVHVQKVTFMLCVGHLEGGKDTCVGDSGGPLMCDGVLQGVTSWGYVP CGTPNKPSVAVRVLVSYVKWIEDTIAENSVDHHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 2mM CaCl ₂ , pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only