

Recombinant Human LAMP2,CD107b (C-6His)

Catalog No: #AP75185



Package Size: #AP75185-1 10ug #AP75185-2 50ug #AP75185-3 500ug #AP75185-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human LAMP2,CD107b (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P13473
Uniprot	P13473
GeneID	3920;
Calculated MW	33.94kDa
Target Sequence	LELNLTDSENATCLYAKWQMNFTVRYETTNTKYKTVTISDHGTVTYNGSICGDDQNGPKIAVQFGPGFSWIAN FTKAASTYSIDSVSFSYNTGDNTTFPDAEDKGILTVDELLAIRIPLNDLFRCNSLSTLEKNDVVQHYWDVLVQAF VQNGTVSTNEFLCDKDKTSTVAPTIHTTVPSPTTTPPKKEPEAGTYSVNNGNDTCLLATMGLQLNITQDKVAS VININPNTTHSTGSCRSHALLRLNSSTIKYLDVFVAVKNENRFYLKEVNISMVYLVNGSVFSIANNLSYWDAPL GSSYMCNKEQTVSVSGAFQINTFDLRVQPFNVTQGGKYSTAQECSLDDDTIVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only