

Brk (Phospho-Tyr447) Polyclonal Antibody

Catalog No: #12268

Package Size: #12268-1 50ul #12268-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Brk (Phospho-Tyr447) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu
Specificity	Phospho-Brk (Y447) Polyclonal Antibody detects endogenous levels of Brk protein only when phosphorylated at Y447.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Brk around the phosphorylation site of Y447.
Target Name	Brk
Modification	Phospho
Other Names	PTK6; BRK; Protein-tyrosine kinase 6; Breast tumor kinase; Tyrosine-protein kinase BRK
Accession No.	Swiss-Prot: Q13882NCBI Gene ID: 5753
Uniprot	Q13882
GeneID	5753;
Target Species	human
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1/500 - 1/2000


Immunofluorescence: 1/200 - 1/1000

ELISA: 1/5000

Not yet tested in other applications

Images

Western Blot analysis of 293 cells using Phospho-Brk (Y447)
Polyclonal Antibody



Note: This product is for in vitro research use only