

Recombinant Human Mannose-Binding Protein C, MBL-2, MBP- C(C-6His)



Catalog No: #AP75189

Orders: order@signalwayantibody.com

Package Size: #AP75189-1 10ug #AP75189-2 50ug #AP75189-3 500ug #AP75189-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Mannose-Binding Protein C, MBL-2, MBP- C(C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P11226
Uniprot	P11226
GeneID	4153;
Calculated MW	25.06kDa
Target Sequence	ETVTCEDAQKTCPAVIACSSSPGINGFPGKDGRDGTKGEKGEPEGQLRGLQGPPGKLGPPGNPGPSGSPGPK GQKGDPGKSPDGDSSLAASERKALQTEMARIKKWLTFSLGKQVGNKFFLTNGEIMTFEKVKALCVKFQASVA TPRNAAENGAIQNLIKEEAFLGITDEKTEGQFVDLTGNRLTYTNWNEGEPNNAAGSDEDCVLLKNGQWNDVP CSTSHLAVCEFPVIDHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 5% Threhalose, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only