

Recombinant Human β -Galactoside α -2,6-Sialyltransferase 1,ST6GAL1 (C-6His)

Catalog No: #AP75201

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75201-1 10ug #AP75201-2 50ug #AP75201-3 500ug #AP75201-4 1mg

Description

Product Name	Recombinant Human β -Galactoside α -2,6-Sialyltransferase 1,ST6GAL1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P15907
Uniprot	P15907
GeneID	6480;
Calculated MW	44.6kDa
Target Sequence	KEKKKGSYYDSFKLQTKFQVLKSLGKLAMGSDSQSVSSSTQDPHRGRQTLGSLRGLAKAKPEASFQVWN KDSSSKNLIPRLQKIWKNYLSMNKYKVSYKGGPGIKFSAEALRCHLRDHVNVSMVEVTDFFPNTSEWEGYLP KESIRTKAGPWGRCAVVSSAGSLKSSQLGREIDDHDAVLRFNAPTANFQQDVGTKTTIRLMNSQLVTTEKRF LKDSLYNEGILIVWDPSVYHSDIPKWYQNPDYNNFFNNTYRKLHPNQPFYILKQMPWELWDILQEISPEEIQ PNPPSSGMLGIIMMMLCDQVDIYEFPLPSKRKTDVCYYYQKFFDSACTMGAYHPLLYEKNLVKHLNQGTDEDIY LLGKATLPGFRTIHCVDHHHHHH
Formulation	Supplied as a 0.2 μ m filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only