

# Recombinant Human Sulfatase Modifying Factor 1,SUMF1 (C-6His)

Catalog No: #AP75202

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75202-1 10ug #AP75202-2 50ug #AP75202-3 500ug #AP75202-4 1mg

## Description

Product Name	Recombinant Human Sulfatase Modifying Factor 1,SUMF1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q8NBK3
Uniprot	Q8NBK3
GeneID	285362;
Calculated MW	38.27kDa
Target Sequence	SQEAGTGAGAGSLAGSCGCGTPQRPGAHGSSAAAHRYRSREANAPGPVPGERQLAHSKMVPIAGVFTMGT DDPQIKQDGEAPARRVTIDAFYMDAYEVSNTFEFEKFNSTGYLTEAEKFGDSFVFEGMLSEQVKTNIQQAVAA APWWLPVKGANWRHPEGPDSTILHRPDHPVLHVSWNDAVAYCTWAGKRLPTEAEWEYSCRGGLHNRLFPW GNKLQPKGQHYANIWQGEFPVTNTGEDGFQGTAPVDAFPPNGYGLYNIVGNAWEWTSDDWWTVHHSVEETL NPKGPPSGKDRVKKGGSYMCHRSYCYRYRCAARSQNTPDSSASNLGFRCAADRLPTMDVDHHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 2mM CaCl <sub>2</sub> , 10% Glycerol, pH 7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only