

Recombinant Human Tryptase  $\beta$ -2,TPSB2 (C-6His)

Catalog No: #AP75205



Package Size: #AP75205-1 10ug #AP75205-2 50ug #AP75205-3 500ug #AP75205-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Tryptase $\beta$ -2,TPSB2 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P20231
Uniprot	P20231
GeneID	64499;7177;
Calculated MW	29.64kDa
Target Sequence	APAPGQALQRVGIVGGQEAPRSKWPWQVSLRVHGPYWMHFCCGSLIHPQWVLTAAHCVGPDVKDLAALRV QLREQHLYYQDQLLPVSRIIVHPQFYTAQIGADIALLEELPEVKVSSHVHTVTLPPASETFPFGMPCWVTGWG DVDNDERLPPPFPLKQVKVPIHENHICDAKYHLGAYTGDDVRIVRDDMLCAGNTRRDSCQGDSSGGPLVCKVN GTWLQAGVVSWECECAQPNRPGIYTRVTYYLDWIHHYVPKPKVDHHHHHH
Formulation	Supplied as a 0.2 $\mu$ m filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only