

Recombinant Human β -2-Microglobulin,B2M (C-6His)

Catalog No: #AP75211



Package Size: #AP75211-1 10ug #AP75211-2 50ug #AP75211-3 500ug #AP75211-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human β -2-Microglobulin,B2M (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P61769
Uniprot	P61769
GeneID	567;
Calculated MW	12.77kDa
Target Sequence	IQRTPKIQVYSRHPAENGKSNFLNCYVSGFHPSDIEVDLLKNGERIEKVEHSDLSFSKDWSFYLLLYTEFTPTE KDEYACRVNHVTLSPKIVKWDRDMVDHHHHHH
Formulation	Lyophilized from a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only