

Recombinant Human Butyrophilin 1A1,BTN1A1 (C-6His)

Catalog No: #AP75213

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75213-1 10ug #AP75213-2 50ug #AP75213-3 500ug #AP75213-4 1mg

Description

| | |
|-----------------|--|
| Product Name | Recombinant Human Butyrophilin 1A1,BTN1A1 (C-6His) |
| Host Species | Human Cells |
| Purification | Greater than 95% as determined by reducing SDS-PAGE. |
| Specificity | Human |
| Accession No. | Q13410 |
| Uniprot | Q13410 |
| GeneID | 696; |
| Calculated MW | 24.98kDa |
| Target Sequence | APFDVIGPPEPILAVVGEDAELPCRLSPNASAEHLELRWFRKKVSPAVLVHRDGREQEAEQMPEYRGRATLV QDGIAGRVALRIRGVRVSDDGEYTCFFREDGSYEEALVHLKVAALGSDPHISMQVQENGEICLECTSVGWY PEPQVQWRTSKGEKFPSTSESERNPDEEGLFTVAASVIIRDTSAKNVSCYIQNLLLQGEKKVEISIPASSLPRVD HHHHHH |
| Formulation | Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2. |
| Storage | Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months. |

Note: This product is for in vitro research use only