

# Recombinant Human C-C Motif Chemokine 3-Like 1,CCL3L1 (C-6His)



Catalog No: #AP75214

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75214-1 10ug #AP75214-2 50ug #AP75214-3 500ug #AP75214-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human C-C Motif Chemokine 3-Like 1,CCL3L1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P16619
Uniprot	P16619
GeneID	414062;6349;
Calculated MW	8.82kDa
Target Sequence	APLAADTPTACCFSYTSRQIPQNFADYFETSSQCSKPSVIFLTRGRQVCADPSEEWVQKYVSDLEPSAVDH HHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only