

# Recombinant Human Ephrin A Receptor 7,EphA7 (C-6His)

Catalog No: #AP75218

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75218-1 10ug #AP75218-2 50ug #AP75218-3 500ug #AP75218-4 1mg

## Description

Product Name	Recombinant Human Ephrin A Receptor 7,EphA7 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q15375
Uniprot	Q15375
GeneID	2045;
Calculated MW	60.19kDa
Target Sequence	QAAKEVLLLD SKAQQTELEWISSPPNGWEEISGLDENYTPIRTYQVCQVMEPNQNNWLR TNWISKGNAQRIF VELKFTLRDCNSLPGVLGTCKETFNLYYETDYDTGRNIRENLYVKIDTIAADESFTQGD LGERKMKLNTEVREI GPLSKKGFYLA FQDVGACIALVSVKVYYKCCWSIIENLAIFPDTVTGSEFSSLVEVRGTCVSSAEEEEAENAPRM HCSAEGEWLVPIGKCICKAGYQQKGDTCPCGRGFYKSSSQDLQCSRCPTH SFSKDEGSSRCECEDGGYRA PSDPPYVACTRPPSAPQNLIFNINQTTVSLEWSPPADNGGRNDVTYRILCKRCSWEQGE CVPCGSNIGYMPQ QTGLEDNYVTMDLLAHANYTFEVEAVNGVSDLSRSQRLFAAVSIT TGQAAPSQVSGVMKERV LQRSVELSW QEPEHPNGVITEYEIKYYEKDQRERTYSTVKTSTASINNLKPGTVVYV FQIRAF TAAGYGNYSRPLDVATLEE ATGKMFEATAVSSEQNPVIVDHHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only