

Recombinant Human Protein Disulfide-Isomerase-Like Protein of the Testis, PDILT (C-6His)



Catalog No: #AP75225

Orders: order@signalwayantibody.com

Package Size: #AP75225-1 10ug #AP75225-2 50ug #AP75225-3 500ug #AP75225-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Protein Disulfide-Isomerase-Like Protein of the Testis, PDILT (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q8N807
Uniprot	Q8N807
GeneID	204474;
Calculated MW	65.54kDa
Target Sequence	SPEVNAGVSSIHITKPVHILEERSLLVLTAPAGLTQMLNQTRFLMVLFHNPSKQSRNLAEEELGKAVEIMGKGN GIGFGKVDITIEKELQQEFGITKAPELKLFFEGNRSEPI SCKGVVESAALVWLRQISQKAF LFNSSSEQVAEFVI SRPLVIVGFFQDLEEEVAELFYDVIKDFPELTFGVITIGNVIGRFHVTLD SVLVFKK GKIVNRQKLINDSTNKQELN RVIKQHLTDFVIEYNTENKDLISELHIMSHMLLFVSKSSESYGIIIQH YKLASKEFQNKILFILVDADEPRNGRVFK YFRVTEVDIPSVQILNLSSDARYKMPSDDITYESLKKFGRSFLSKNATKHQSSEEIPKYWDQGLVKQLVGKNFN VVVFDKEKDVFMFYAPWSKCKMLFPLLEELGRKYQNHSTIIIAKIDVTANDIQLMYLDRYPFFRLFPSGSQQ AVLYKGEHTLKGFSDFLESHIKTKIEDEDELLSVEQNEVIEEEVLAEEKEVPMMRKGLPEQQSPELENMTKYVS KLEEPAGKKKTSEEVVVVAKPKGPPVQKKKPKVKEELVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 10% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only