

# Recombinant Human Nogo-66 Receptor-Related 3,NgR3,RTN4RL1 (C-6His)



Catalog No: #AP75229

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75229-1 10ug #AP75229-2 50ug #AP75229-3 500ug #AP75229-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Nogo-66 Receptor-Related 3,NgR3,RTN4RL1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q86UN2
Uniprot	Q86UN2
GeneID	146760;
Calculated MW	45.49kDa
Target Sequence	CPRDCVCYPAPMTVSCQAHNFAAIEGIPVDSERVFLQNNRIGLLQPGHFSPAMVTLWIYSNNITYIHPSTFEG FVHLEELDLGDNRQLRTLAPETFQGLVKLHALYLYKCGLSALPAGVFGGLHSLQYLYLQDNHIEYLQDDIFVDL VNLSHLFLHGNKLWSLPGTFRGLVNLDRLLLHENQLQWVHHKAFHDLRRLTTLFLFNNSLSELGCECLAPLG ALEFLRLNGNPWDCGCRARSLWEWLQRFRRGSSSAVPCVSPGLRHGQDLKLLRAEDFRNCTGPASPHQIKSH TLTTTDRAARKEHHSHPGPTRSKGPHGPRPGHRKPGKNCTNPRNRNQISKAGAGKQAPELPDYAPDYQHK FSFDIMPTARPKRKGKCARRTPIRAPSGVQQAVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 5% Threhalose, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only