

Recombinant Human Anterior Gradient Protein 3 Homolog,AG-3,BCMP11,AGR3 (C-6His)

Catalog No: #AP75249

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75249-1 10ug #AP75249-2 50ug #AP75249-3 500ug #AP75249-4 1mg

Description

Product Name	Recombinant Human Anterior Gradient Protein 3 Homolog,AG-3,BCMP11,AGR3 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q8TD06
Uniprot	Q8TD06
GeneID	155465;
Calculated MW	18.04kDa
Target Sequence	IAIKKEKRPPQTLSRGWGDDITWVQTYEEGLFYAQKSKKPLMVIHHLEDCQYSQALKKVFAQNEEIQEMAQNK FIMLNLMHETTDKNLSPDGQYVPRIMFVDPSTVRADIAGRYSNRLYTYEPRDLPLLIENMKKALRLIQSELVDH HHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 2mM EDTA, pH 8.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only