

Recombinant Human Alkaline Phosphatase,ALPL (C-6His)



Catalog No: #AP75250

Orders: order@signalwayantibody.com

Package Size: #AP75250-1 10ug #AP75250-2 50ug #AP75250-3 500ug #AP75250-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Alkaline Phosphatase,ALPL (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P05186
Uniprot	P05186
GeneID	249;
Calculated MW	54.5kDa
Target Sequence	LVPEKEKDPKYWRDQAQETLKYALELQKLNTNVAKNVIMFLGDGMGVSTVTAARILKQQLHHNPGEETRLEM DKFPFVALSKTYNTNAQVPDSAGTATAYLCGVKANEGTVGVSAATERSRCNTTQGNEVTSILRWAKDAGKSV GIVTTTRVNHATPSAAYAHSADRDWYSDNEMPPEALSQGCKDIAYQLMHNIRDIDVIMGGGRKMYMPKNKTD VEYESDEKARGTRLDGLDLVDTWKSFKPRYKSHSHFIWNRTELLTLDPHNVDYLLGLFEPGDMQYELNRNNVT DPSLSEMVVVAIQILRKNPKGFLLVEGGRIDHGHHEGKAKQALHEAVEMDRAIGQAGSLTSSSEDTLTVVTAD HSHVFTFGGYTPRGNSIFGLAPMLSDTDKPKFTAILYGNGPGYKVVGGGERENVSMVDYAHNNYQAQSAVPLR HETHGGEDVAVFSKGPMALLHGVHEQNYVPHVMAYAACIGANLGHCAPASSVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM HEPES, 150mM NaCl, 2mM MgSO ₄ , 0.1mM ZnCl ₂ , pH 7.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only