

# Recombinant Human Butyrylcholine Esterase, BCHE (C-6His)



Catalog No: #AP75254

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75254-1 10ug #AP75254-2 50ug #AP75254-3 500ug #AP75254-4 1mg

## Description

Product Name	Recombinant Human Butyrylcholine Esterase, BCHE (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P06276
Uniprot	P06276
GeneID	590;
Calculated MW	66.12kDa
Target Sequence	EDDIIIATKNGKVRGMNLTVFGGTVTAFLGIPYAQPPLGRLRFKPKQSLTKWSDIWNATKYANSCCQIDQSFP GFHGSEMWNPNLTDLSEDCLYLNVWIPAPKPKNATVLIWIYGGGFQTGTSSLHVYDGKFLARVERVIVSMNY RVGALGFLALPGNPEAPGNMGLFDQQLALQWVQKNIAAFGGNPKSVTLFGESAGAASVSLHLLSPGSHSLFT RAILQSGSFNAPWAVTSLYEARNRNLNLAKLTGCSRENETEIIKCLRNDKDPQEILLNEAFVVPYGTPLSVNFGPT VDGDFLTDMPDILLELGQFKKTQILVGVNKDEGTAFLVYGAPGFSKDNNSIITRKEFQEGLKIFFPGVSEFGKES ILFHYTDWVDDQRPENYREALGDVVDYDFICPALEFTKKFSEWGNNAFFYYFEHRSSKLPWPEWMGMHG YEIEFVFGLEPLERRDNYTKAEELSRIVKRWANFAKYGNPNETQNNSTSWPVFKSTEQKYLTLNTESTRIMTK LRAQQCRFWTSFFPKVLEMTGNIDEAEWEWKAGFHRWNNYMMDWKNQFNDYTSKKESCGLVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only