

## Recombinant Human CD19 (C-Fc)

Catalog No: #AP75257



Package Size: #AP75257-1 10ug #AP75257-2 50ug #AP75257-3 500ug #AP75257-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CD19 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P15391
Uniprot	P15391
GeneID	930;
Calculated MW	57.3kDa
Target Sequence	PEEPLVVKVEEGDNAVLQCLKGTS DGPTQQLTWSRESPLKPFLLKLSLGLPGLGIHMRPLAIWLFIFNVSQQMG GFYLCQPGPPSEKAWQPGWTVNVEGSGELFRWNVSDLGGLGCGLKNRSSEGPSSPSGKLMSPKLYVWAK DRPEIWE GEPPLPPRDSL NQSLSQDLTMAPGSTLWLS CGVPPDSVSRGPLSWTHVHPKGP KSLLSLELKDD RPARDMWMETGLLLPRATAQDAGKYYCHRG NLTMSFHLEITARPVLWHWLLRTGGWKVDDIEGRMDEPKS CDKTHTCP PCAPPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKP REEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQV SLTCLVKG FYPSPDI AVEWESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNV FSCSVMHEALHNH YTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only