

Recombinant Human CD3  $\epsilon$ ,CD3E (C-6His)

Catalog No: #AP75260



Package Size: #AP75260-1 10ug #AP75260-2 50ug #AP75260-3 500ug #AP75260-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CD3 $\epsilon$ ,CD3E (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P07766
Uniprot	P07766
GeneID	916;
Calculated MW	12.79kDa
Target Sequence	DGNEEMGGITQTPYKVISISGTTVILTCPQYPGSEILWQHNDKNIGGDEDDKNIGSDEDHLSLKEFSELEQSGYY VCYPRGSKPEDANFYLYLRARVCENCMEMDVDHHHHHH
Formulation	Lyophilized from a 0.2 $\mu$ m filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only