

IL-13R α 1 (Phospho-Tyr405) Polyclonal Antibody

Catalog No: #12276

Package Size: #12276-1 50ul #12276-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-13R α 1 (Phospho-Tyr405) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	Phospho-IL-13R α 1 (Y405) Polyclonal Antibody detects endogenous levels of IL-13R α 1 protein only when phosphorylated at Y405.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human IL-13R α 1 around the phosphorylation site of Y405.
Target Name	IL-13R α 1
Modification	Phospho
Other Names	IL13RA1; IL13R; IL13RA; Interleukin-13 receptor subunit alpha-1; IL-13 receptor subunit alpha-1; IL-13R subunit alpha-1; IL-13R-alpha-1; IL-13RA1; Cancer/testis antigen 19; CT19; CD antigen CD213a1
Accession No.	Swiss-Prot: P78552NCBI Gene ID: 3597
Uniprot	P78552
GeneID	3597;
Target Species	human
SDS-PAGE MW	49kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1/500 - 1/2000


Immunofluorescence: 1/200 - 1/1000

ELISA: 1/10000

Not yet tested in other applications

Images

Western Blot analysis of 3T3 cells using Phospho-IL-13R α 1 (Y405) Polyclonal Antibody



Note: This product is for in vitro research use only