

Recombinant Human Carboxypeptidase A1, CPA1 (C-6His)

Catalog No: #AP75273

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75273-1 10ug #AP75273-2 50ug #AP75273-3 500ug #AP75273-4 1mg

Description

Product Name	Recombinant Human Carboxypeptidase A1, CPA1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P15085
Uniprot	P15085
GeneID	1357;
Calculated MW	46.57kDa
Target Sequence	KEDFVGHQVLRISVADEAQVQKVKELLEDLEHLQLDFWRGPAHPGSPIDVRVPFPSIQAVKIFLESHGISYETMIE DVQSLLDDEEQEQMFAFRSRARSTDFTFNATYHTLEEIYDFLDLLVAENPHLVSKIQIGNTYEGRPIYVLKfstgg SKRPAIWIDTGIHSREWVTQASGVWFakkITQDYGDAAFTAILDTLDIFLEIVTNPdGFafthSTNRMWRKTR SHTAGSLCIGVDPNRNWDAGFGLSGASSNPCSETYRgKfANSEVEVKSIVDFVKDHGNIKAFISIHSYSQLLMY PYGYKTEPVPDQDELQLSKAAVTALASLYGtKfNYGSIKAIYQASGSTIDWTYSQGIKYSFTfELRDTGRYGF LLPASQIIPTAKETWLALLTIMEHTLNHPYVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mm NaCl, pH 7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only