

Recombinant Human Carboxypeptidase A2, CPA2 (C-6His)

Catalog No: #AP75274

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75274-1 10ug #AP75274-2 50ug #AP75274-3 500ug #AP75274-4 1mg

Description

Product Name	Recombinant Human Carboxypeptidase A2, CPA2 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P48052
Uniprot	P48052
GeneID	1358;
Calculated MW	45.91kDa
Target Sequence	LETFVGDQVLEIVPSNEEQIKNLLQLEAQEHLQLDFWKSPPTTPGETAHVRVPFVNVQAVKVFLGSQGIAYSIMIE DVQVLLDKENEEMLFNRRRERSGNFNFGAYHTLEEISQEMDNLVAEHPGLVSKVNIGSSFENRPMNVLKFST GGDKPAIWL DAGIHAREWVTQATALWTANKIVSDYGKDPSITSILDALDIFLLPVTNPDGYVFSQTKNRMWRKT RSKVSGSLCVGVDPNRNWDAGFGGPGASSNPCSDSYHGPSANSEVEVKSIVDFIKSHGKVKAFITLHSYSQL LMFPYGYKCTKLDDFDELSEVAQKAAQSLRSLHGTYKVGPICSVIYQASGGSIDWSYDYGIKYSFAFELRDT GRYGFLPARQILPTAETWLGLKAIMEHVRDHPYVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 150mm NaCl, pH 7.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only