

Recombinant Human Carboxypeptidase A4, CPA4 (C-6His)

Catalog No: #AP75275

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75275-1 10ug #AP75275-2 50ug #AP75275-3 500ug #AP75275-4 1mg

Description

Product Name	Recombinant Human Carboxypeptidase A4, CPA4 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9UI42
Uniprot	Q9UI42
GeneID	51200;
Calculated MW	46.63kDa
Target Sequence	GQEKFFGDQVLRINVRNGDEISKLSQLVNSNNLKLNFWKSPSSFNRPVDVLVPSVSLQAFKSFLRSQGLEYAV TIEDLQALLDNEDEMQRHNEGQERSSNNFNYGAYHSLEAIYHEMDNIAADFPDLARRVKIGHSFENRPMYVLK FSTGKGVRRPAVWLNAGIHSREWISQATAIWTARKIVSDYQRDPAITSILEKMDIFLLPVANPDGYVYTQTQNR LWRKTRSRNPGSSCIGADPNRNWNASFAGKGGASDNPCSEVYHGHANSEVEVKSVDVFIQKHGNFKGFIDL HSYSQLLMYPYGYSVKKAPDAEELDKVARLAALKALASVSGTEYQVGPTCTTVYPASGSSIDWAYDNGIKFAFT FELRDTGTYGFLLPANQIIPTAEETWLGLKTIMEHVRDNLVVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only