

# Recombinant Human Corticotropin-Releasing Factor-Binding Protein, CRHBP (C-6His)



Catalog No: #AP75276

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75276-1 10ug #AP75276-2 50ug #AP75276-3 500ug #AP75276-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Corticotropin-Releasing Factor-Binding Protein, CRHBP (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P24387
Uniprot	P24387
GeneID	1393;
Calculated MW	34.38kDa
Target Sequence	YLELREAADYDPFLLFSANLKRELAGEQPYRRALRCLDMLSLQGQFTFTADRPQLHCAAFFISEPEEFITIH YDQVSDICQGGDFLKVFDGWILKGEKFPSSQDHPLPSAERYIDFCESGLSRRSIRSSQNVA MIFFRVHEPGNGFTLTIKTDPNLFPCNVISQTPNGKFTLVVPHQHRNCSFSIIPVVIKISDLTLGHV NGLQLKKSSAGCEGIGDFVELLGGTGLDPSKMTPLADLCYPFHGPAQMKVGC DNTVVRMVSSGKHVNRVTFEYRQLEPYELENPNGNSIGEFCL SGLVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 150mm NaCl, pH 7.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only